

Setaflash™ Series 3

The safest choice

Approved in over 1000 international product specifications and regulations

- ISO 3679; ISO 3680 (obs); ISO 9038
- ASTM D3278; D3828; D4206; D7236; D8174; EPA 1020 B; SW-846
- IP 523; 524; 534; 602; BS 2000-523



1 minute test • 2ml sample • 0 to 300°C • Portable • Conforms



Why use Setaflash Series 3?

- Conforms to international regulations
- 1 or 2 minutes test time
- 2ml sample
- 0°C to 300°C

Key features

- Compact, portable and rugged design
- Best published precision of any flash point method
- Suitable for unknown samples using ramp mode
- Automated flash detection
- Automatic barometric pressure correction
- USB port and result storage
- Low cost of purchase and maintenance



Where can I find more information?

Learn more about Setaflash testing:

www.stanhope-seta.co.uk/small-scale-flashpoint-testing.asp



Watch the Series 3 video demonstration:

www.stanhope-seta.co.uk/5134/Setaflash-Series-3-Flash-Point-Tester




Operator sequence

4 easy steps:



Select test temperature



Inject 2ml sample



Press 




Dip the test flame, flash detection is automatic



Trusted Precision and Accuracy

- The most precise flash point test
- Proven performance
- Approved in an extensive range of specifications and regulations



USB Data Transfer

- 1GB data capacity, holding 100,000 results
- Results can be saved to a USB stick



Calibration and Verification

- Software enables the user to calibrate temperature and barometric pressure
- Verification through Seta's Small Scale Certified Flash Point material (reference 99878-3)



Small Sample - 2ml

- Small sample size of only 2ml
- Reduces costs and waste per test



Service and Maintenance

- Rugged design
- Minimum maintenance
- Visit our technical support website: www.stanhope-seta.co.uk/support



Rapid Results in just 1 or 2 Minutes

- Fastest method of checking flash point
- Results determined in just one or two minutes
- Maximises operability



Search for the Flash Point

- Test unknown samples using ramp mode
- Ramp mode is integral to 33000-2, 33200-3 and 33250-3 models
- User adjustable ramp rates dip interval

| Technical Specifications: | | | |
|----------------------------|---|--|--|
| Seta Part Number: | 30000-2 | 33000-2 | 33200-3 & 33250-3 |
| Cup: | Closed | Closed | Closed |
| Material: | Aluminium | Aluminium | Aluminium or stainless steel insert |
| Temperature Range: | Ambient to 300°C | Ambient to 300°C | 10 to 135°C |
| Temperature Ramp Mode: | N/A | Automatic or custom settings (ramp rate: up to 6°C/min and beep rate: up to 5°C/B intervals) Extended test range, 100°C from start temperature | Automatic or custom settings (ramp rate: up to 6°C/min and beep rate: up to 5°C/B intervals) Extended test range, 100°C from start temperature |
| Temperature Selection: | °C or °F | °C or °F | °C or °F |
| User Interface: | Colour LCD display - touch screen sensitive | Colour LCD display - touch screen sensitive | Colour LCD display - touch screen sensitive |
| Navigation: | Control knob (rotational and push-down) | Control knob (rotational and push-down) | Control knob (rotational and push-down) |
| Data Storage and Download: | 1GB (internal memory and integrated USB port) | 1GB (internal memory and integrated USB port) | 1GB (internal memory and integrated USB port) |
| Fuel Supply: | Integrated gas tank (removable) | Integrated gas tank (removable) | Integrated gas tank (removable) |
| Sample Size: | 2ml or 4ml | 2ml or 4ml | 2ml or 4ml |
| Pressure Correction: | Automatic barometric correction | Automatic barometric correction | Automatic barometric correction |
| Voltage Supply: | 110 to 250V a.c.-50/60 Hz (universal) | 110 to 250V a.c.-50/60 Hz (universal) | 110 to 250V a.c.-50/60 Hz (universal) |
| Power: | 200W (Max.) | 200W (Max.) | 200W (Max.) |

Products & Specifications

Diesel | Aviation Fuel | Marine Fuel | Biodiesel | FAME | Chemicals | Pharmaceuticals | Paints | Cosmetics | Waste Flavours | Inks | Waxes | Adhesives | Oils - Lubricating, Hydraulic, Base, Mineral, Used/Cooking Oils | Pastes

Applications

| | |
|-----------------|---------------------------------|
| QC | |
| Quality Control | Transport & Storage Regulations |
| | |
| CLP Regulations | Waste Disposal Regulations |

| Application: | Test Method: | Who Says So: |
|-----------------------|---------------------------------|--|
| Transport Regulation | Small Scale, other Closed Cups | UN GHS; IATA; ADR; IMDG; CLP; DOT CFR 49-173.120 |
| Aviation Turbine Fuel | Small Scale, Abel, Tag | ASTM D1655; Def Stan 91-91; IATA Guidance Material; ATA 103; AFQRJOS |
| Gas Turbine Fuel | Small Scale, Pensky-Martens | ASTM 2880 |
| Diesel Fuel | Small Scale, Pensky-Martens | ASTM D975; D7467 |
| Kerosines | Small Scale, Tag | ASTM D3699 |
| Biodiesel | Small Scale | EN14213; EN14214; ASTM D6751 |
| Fuel Oil | Small Scale, Pensky-Martens A&B | ASTM D396; ISO 8217 |
| Naphthas | Small Scale, Tag | ASTM D3734; D3735 |
| Health & Safety | Small Scale | OSHA 29 CFR 1910.106 & 1200; CPSC CFR 16-1500 43a |
| Water Borne Paints | Small Scale | ISO 3679; ISO 3680 |
| Waste Products | Small Scale | EPA 1020 B; CFR40 261.21; ASTM D8174; 2008/98/EC |

“For us the Setaflash Tester is a valuable piece of kit and we’d struggle without it. Our main use is to detect whether a waste is hazardous or not, as this can have knock on effects in terms of cost, safety, compliance and potential processing routes. It’s very useful to have it at hand as sending samples out for off site analysis simply would not be an option at times, due to potentially lengthy turn around times.”

*Nick Richardson,
Laboratory Supervisor,
Castle Environmental*

Discover the entire Setaflash Range

| Seta Model: | Seta Part No: | Temperature Range: | Heating/Cooling Method: |
|---|---------------|---|---|
| Series 3 Activecool | 33200-3 | 10 to 135°C | Peltier Cell |
| Series 3 Activecool (Corrosion Resistant) | 33250-3 | 10 to 135°C | Peltier Cell |
| Series 3 'Plus' (Auto Ramp) with touchscreen | 33000-2 | Ambient to 300°C *(0 to 300°C with coolant module) | Cartridge Heater |
| Series 3 with touchscreen | 30000-2 | Ambient to 300°C *(0 to 300°C with coolant module) | Cartridge Heater |
| Series 3 Open Cup | 31000-0 | Ambient to 300°C *(0 to 300°C with coolant module) | Cartridge Heater |
| Series 8 High Temperature (Hot Wire) | 82000-0 | Ambient +5 to 300°C | Ceramic Pad, Forced Air (post-test cool down) |
| Series 8 High Temperature (Gas) | 82050-0 | Ambient +5 to 300°C | Ceramic Pad, Forced Air (post-test cool down) |
| Series 8 Activecool (Hot Wire) | 82100-0 | Air: 10 to 130°C Water assisted: -20 to 130°C | Peltier cell |
| Series 8 Activecool (Gas) | 82110-0 | Air: 10 to 130°C Water assisted: -20 to 130°C | Peltier cell |
| Series 8 Activecool (Corrosion Resistant, Hot Wire) | 82150-0 | Air: 10 to 130°C Water assisted: -20 to 130°C | Peltier cell |
| Series 8 Activecool (Corrosion Resistant, Gas) | 82160-0 | Air: 10 to 130°C Water assisted: -20 to 130°C | Peltier cell |



For more information please visit: www.stanhope-seta.co.uk