

1



1 LLG-uniTOPDISPENS, Dispensers for bottles

With glass piston and clear glass cylinder.
 Accurate and reliable volume control
 Visible dosing procedure
 PP valves block providing optimized protection of the glass dosing cylinder
 High precision of the glass piston to have a perfect seal in the glass cylinder
 Ergonomic fingers-protector
 Supplied with 2 adapters in PP: GL40 and GL45
 Accuracy Max. Vol. 0.6 ($\leq \pm A\%$)
 Precision Max. Vol. 0.2 ($\leq CV\%$)

Type	Sample grad. volume mL mL	PK	Cat. No.
LLG-uniTOPDISPENS	2 - 10 0,2	1	9283806
LLG-uniTOPDISPENS	10 - 50 1,0	1	9283808

2



9055513

2 Syringe filters & filter paper

Description	PK	Cat. No.
SYRINGE FILTER Ø 25 mm, CA, 0,45 µm, yellow. STERILE	50	9055513
SYRINGE FILTER Ø 25 mm, CA, 0,45 µm, yellow	500	9055503
SYRINGE FILTER Ø 25 mm, Nylon, 0,45 µm, transparent	500	9055523
SYRINGE FILTER Ø 25 mm, Nylon, 0,45 µm, transparent	50	6272641
SYRINGE FILTER Ø 25 mm, Nylon, 0,45 µm, transparent	45	6272643
FILTER PAPER, plan disc, Ø 150 mm RAPID	100	4006418
FILTER PAPER, plan disc, Ø 185 mm RAPID	100	4006419
FILTER PAPER, plan disc, Ø 210 mm RAPID	100	3997503
FILTER PAPER, plan disc, Ø 250 mm RAPID	100	4006420
FILTER PAPER, folded disc, Ø 150 mm EXTRA-RAPID	100	4006219
FILTER PAPER, folded disc, Ø 190 mm EXTRA-RAPID	100	4006220
FILTER PAPER, folded disc, Ø 210 mm EXTRA-RAPID	100	4006221
FILTER PAPER, folded disc, Ø 250 mm EXTRA-RAPID	100	4006223
FILTER PAPER, in ream, 500x500 mm BIBULA, 67 g for rough impurities (sheets)	100	G01526
FILTER PAPER, in ream, 500x500 mm, 67 g for rough impurities (sheets)	500	G01989

3



3M respirator

3 DPI/PPE - Personal Protection Devices

Description	PK	Cat. No.
Protection spectacle "Classic" blue frame, clear lenses	1	9006001
Respirator, Type 8710, FFP 1 without valve	20	6050640
Respirator Type 8810, FFP 2 without valve	20	6078942
Respirator Type 8812, FFP 1 with valve	10	6202080
Respirator Type 8822, FFP 2 with valve	10	6052910
Disposable gloves, Standard, Nitrile, size S, powder free	100	9006375
Disposable gloves, Standard, Nitrile, size M, powder free	100	9006376
Disposable gloves, Standard, Nitrile, size L, powder free	100	9006377
Disposable gloves, Standard, Nitrile, size XL, powder free	100	9006378
Protective coat, woman, model Carla, white, elastic wrists, size M	1	YL0600201
Protective coat, woman, model Carla, white, elastic wrists, size L	1	YL0600202
Protective coat, man, model Capri, white, elastic wrists, size M	1	YL0600211
Protective coat, man, model Capri, white, elastic wrists, size L	1	YL0600212
Protective coat, man, model Capri, white, elastic wrists, size XL	1	YL0600213

4



TFI-260 thermometer

4 Consumables

Description	PK	Cat. No.
Wash bottle with extractable pin, 500 mL	1	YK0000183
Safety pipetting ball, universal type, red	1	6270266
Bunsen Poligas lamp, with thermocouple and tap	1	YL0084000
Absorbent paper coil, double thickness wadding. 625 rips. 300x400 mm	2	YL0325500
Absorbent paper coil, 800 rips. 80% rips, 20% recycled	2	YL0325750
Foldable digital thermometer TLC 700, -30 to 220 °C	1	9908018
Min/Max Thermometer with outdoor sensor.	1	9243165
Inside -20 to 70 °C, outside -50 to 70 °C		
Infrared thermometer TFI-260. Measuring range -60 to 550 °C	1	6272762

1 Glassware

Borosilicate glass 3.3

Description	PK	Cat. No.
Beaker with spout, low form, 100 mL	1	YA0014505
Beaker with spout, low form, 250 mL	1	YA0014705
Beaker with spout, low form, 600 mL	1	YA0014905
Beaker with spout, low form, 800 mL	1	YA0015005
Flask, narrow mouth, 100 mL	1	YA0180405
Flask, narrow mouth, 250 mL	1	YA0180805
Flask, narrow mouth, 500 mL	1	YA0181105
Filtration flask with spout, 1000 mL	1	YA0186102
Graduated cylinder with spout, class A, 50 mL	1	YA0089010
Graduated cylinder with spout, class A, 100 mL	1	YA0089210
Graduated cylinder with spout, class A, 250 mL	1	YA0089410
Volumetric flask, class A, NS 12/21, with Moplen cap, capacity 100 mL	1	YA0210900
Volumetric flask, class A, NS 14/23, with Moplen cap, capacity 200 mL	1	YA0211000
Volumetric flask, class A, NS 19/26, with Moplen cap, capacity 500 mL	1	YA0211200
Funnel, short shaft, diam. 75 mm	1	YA0158205
Funnel, short shaft, diam. 100 mm	1	YA0159002
Volumetric pipette, 1-mark graduation, class A, capacity 5 mL	1	YA0277400
Volumetric pipette, 1-mark graduation, class A, capacity 10 mL	1	YA0277500
Volumetric pipette, 1-mark graduation, class A, capacity 20 mL	1	YA0277700
Volumetric pipette, 1-mark graduation, class A, capacity 25 mL	1	YA0277800
Graduated pipette, class A, capacity 5 mL, div. 1/10	1	YA0257800
Graduated pipette, class A, capacity 10 mL, div. 1/10	1	YA0258000
Mohr burette, Schelbach glass, class A, capacity 25 mL, div. 1/10	1	YA0052400
Pellet burette, Schelbach glass, class A, capacity 25 mL, div. 1/10	1	YA0065300
Pellet burette, amber glass, class A, capacity 25 mL, div. 1/10	1	YA0065400



2 3 Plastic consumables

Description	PK	Cat. No.
Measuring jug, low form, 500 mL	1	YK0001156
Measuring jug, low form, 1000 mL	1	YK0001157
Measuring jug, low form, 2000 mL	1	YK0001158
Measuring jug, low form, 3000 mL	1	YK0001159
Measuring jug, low form, 5000 mL	1	YK0001160
Measuring jug, tall form, 250 mL	1	YK0001086
Measuring jug, tall form, 500 mL	1	YK0001087
Measuring jug, tall form, 1000 mL	1	YK0001088
Measuring jug, tall form, 2000 mL	1	YK0001089
Measuring jug, tall form, 3000 mL	1	YK0001090
Bottle with graduated dispenser. Measuring container 25 mL; bottle 250 mL	1	YK0001475
Bottle with graduated dispenser. Measuring container 25 mL; bottle 500 mL	1	YK0001476
Bottle with graduated dispenser. Measuring container 25 mL; bottle 1000 mL	1	YK0001477
Bottle with graduated dispenser. Measuring container 50 mL; bottle 250 mL	1	YK0001480
Bottle with graduated dispenser. Measuring container 50 mL; bottle 500 mL	1	YK0001481
Bottle with graduated dispenser. Measuring container 50 mL; bottle 1000 mL	1	YK0001482
Graduated cylinder, low form, 10 mL	1	YK0000814
Graduated cylinder, low form, 25 mL	1	YK0000815
Graduated cylinder, low form, 50 mL	1	YK0000816
Graduated cylinder, low form, 100 mL	1	YK0000817
Graduated cylinder, low form, 250 mL	1	YK0000820
Graduated cylinder, low form, 500 mL	1	YK0000821
Graduated cylinder, low form, 1000 mL	1	YK0000822
Graduated cylinder, low form, 2000 mL	1	YK0000823
Funnel, diam. 66 mm	1	YK0000148
Funnel, diam. 81 mm	1	YK0000150
Funnel, diam. 120 mm	1	YK0000153
Vertical pipette stand. 18 small and 10 large holes places	1	YK0000265
Circular base with PVC rod, for art. YK0000265	1	YK0000262
Vertical pipette stand, 200x75 mm	1	YK0000255
Pipette stand, PP, vertical and rotatable, 94 places	1	9273895



Measuring jug, low form



Bottle with dispenser