

1 2 3 4 5 Afrometers

For CO₂ detection.

Description	PK	Cat. No.
PIN AFROMETER, for cork caps, 0-6 bar pressure gauge	1	A113501
PIN AFROMETER with tap, for cork caps, 0-6 bar pressure gauge	1	A113502
PIN VACUUM GAUGE with tap, for cork caps, -1+5 bar pressure gauge	1	A113502/V
"UNIVERSAL" AFROMETER with case, for any cap, 0-6 bar pressure gauge	1	A113503
DIGITAL "UNIVERSAL" AFROMETER with case, for any cap, 0-6 bar pressure gauge	1	A113504
AFROMETER for pressure measurement during the fermentation in bottle	1	A113506



A113503




A113504



A113506



A113501



A113502 - A113502/V

Accessories and spares for afrometers

Description	PK	Cat. No.
Spare pin	1	X200196
Ejecting tip	1	X200197
Pressure gauge 0 - 6 Bar	1	X200194
Pressure gauge - 1 / 5 Bar	1	X200542
Digital pressure gauge 0 - 6.000, Division 0.001 Bar	1	X200781
Pressure gauge 0 - 6 Bar, with S.I.T. Certificate	1	X200194C
Pressure gauge - 1 / 5 Bar, with S.I.T. Certificate	1	X200542C
Silicon seal for Universal Afrometer	1	X200137

6 WITHDRAWING DEVICE and AFROMETER for sparkling wine

Device that allows bottles to be filled with sparkling wine in order to directly sample on the autoclave under pressure.

Description	PK	Cat. No.
WITHDRAWING DEVICE and AFROMETER for sparkling wine	1	A11348

6



7 CARBODOSEUR for CO₂ detection

100 mL cylinder with sealed screw cap and vent hose.

Description	PK	Cat. No.
Carbodoseur	1	A13001

7



8 9 Oxygen meters

Instruments for measuring oxygen dissolved (and other parameters) in wines.

Description	PK	Cat. No.
HQ30d - Digital portable instrument, 1 channel for pH, conductivity, dissolved oxygen & redox. With pH electrode	1	39987198
HQ30d - Digital portable instrument, 1 channel	1	39987101
LDO luminescent probe for standard oxygen analysis, with 3m cable	1	3998763

8



HQ30D

9



Probe